



Submittal Data Sheet

Project

Information

Project Number _____ Approval _____

- **Aging medical gas outlets?**
- **Leakage problems?**
- **Need increased flow?**
- **Can't find parts?**

Convert your facility from many styles of outlets to Tri-Tech Medical's Ohio/Ohmeda Diamond III style outlet.

There is no need to cut existing outlets out of the wall. Leaking, obsolete, or discontinued outlets can be quickly converted at a significantly lower cost than replacing them with new construction outlets.



Frontall™ feature allows service of primary check valve from the front without removing the faceplate !

Specification:

The medical gas outlet conversion kit shall be designed so that, once installed, routine service of the primary check valve can be accomplished **without** removing the nameplate or gas specific portions. So that the primary check valve can be removed for service without shutting off the gas supply to the outlet, the secondary check valve shall operate automatically to stop the free flow of the pressure gas. Medical gas outlets that require the removal of the nameplate or gas specific components for routine service shall not be acceptable. The outlet nameplate shall be permanently color coded and sealed behind an easy to clean plastic protective cover.

Features:

- Accepts Ohmeda/Ohio gas specific compatible male adaptors
- Unique patented front loaded feature
- Ease of maintenance using an insertion tool
- Gas specific indexing
- Modular design
- All outlets are 100% tested for flow & leakage
- Full faceplate color coding
- Made in the U.S.A.

How to Order:

Existing outlet description	N2O	Air	Vac	O2
Melco/Foregger/Linde 9/16-18	Q3104	Q3116	Q3122	Q3124
Oxequip 07, Medstar	Q3704	Q3716	Q3722	Q3724
Chemetron/NCG 5/8-27	Q3804	Q3816	Q3822	Q3824
Chemetron/NCG 3/4-16	Q3904	Q3916	Q3922	Q3924
Ohmeda Diamond I	Q4004	Q4016	Q4022	Q4024
Ohmeda Diamond II/III	Q4104	Q4116	Q4122	Q4124
Ohio "40" Series	Q4204	Q4216	Q4222	Q4224
Puritan "B" series	Q4604	Q4616	Q4622	Q4624
Amico	Q4704	Q4716	Q4722	Q4724
Hill-Rom DiamondCare™	Q4804	Q4816	Q4822	Q4824

- Trim plates are available (at an extra cost) for multiple outlets and can be installed against the wall behind the outlet faceplate. Trim plates are custom made based on the distance between outlet centers.
- Contact Tri-Tech Medical for availability of additional families of conversion kits and/or gas services