



SCS-2000 Suction Roll Stand

The SCS-2000 Model has a ball-pivot holder for attachment of the suction regulator. This allows movement of the regulator in any direction making visual observation of the suction rate easier. The four-port, on-off valve allows quick change from canister-to-canister during the suctioning process. The universal canister holder can accommodate numerous canister configurations.

All Roll Stands are sold without suction regulators, suction tubing/hydrophobic filters, suction canisters, or canister holders.

Ordering information:

60-0601B

Continuous 3-mode Suction Regulator 0-200 mm Hg and Full Vacuum, no connections (1/8" npt female connections)

503477

Suction tubing / hydrophobic in-line filter assembly for SCS2000

64-0602

Plastic Ring for 1200 cc disposable suction canister

64-0612

1200 cc Disposable Suction Canister

64-0606

Stainless Steel mounting bracket for reusable suction bottle and cap assembly

64-0603

½ gallon (1200 cc) Reusable Polycarbonate Suction Bottle and Cap Assembly

SCS2000

Suction Roll Stand – Includes ball-pivot holder, four-port on-off valve, and slide bracket

