



DISS Quick Connects

Function:

DISS quick connect fittings (also known as “Zip Nuts”) are designed to be used anywhere you would normally use a wrench or hand-tight style DISS fitting. To attach you simply push the quick connect onto the outlet – then turn the nut clockwise until it is tight. (Usually only about 1/3rd of a turn). The advantage of the quick connect is that it attaches to the outlet much more quickly than the wrench or hand-tight style DISS nut and nipples. The DISS quick connector also prevents operator fatigue with limited wrist motion that may help prevent carpal tunnel syndrome.

Features:

- Available for all seven common medical gas services
- Fast to connect – saves you time
- O-ring seal
- Color coded and labeled for specific gas services
- Available in NPT and hose barb styles
- Large hand tight size nut for “hands only” Assembly (no wrench required)
- Safe – working pressure rated for over 5 times the normal working pressure



How to Order:

DISS x 1/8 NPT Male

Part #	Gas Service	Description
Z5004-2	Nitrous Oxide	QC/DISS x 1/8 NPT Male
Z5012-2	Nitrogen	QC/DISS x 1/8 NPT Male
Z5016-2	Air	QC/DISS x 1/8 NPT Male
Z5022-2	Vacuum	QC/DISS x 1/8 NPT Male
Z5024-2	Oxygen	QC/DISS x 1/8 NPT Male
Not available	Evacuation	QC/DISS x 1/8 NPT Male

DISS x 1/4 NPT Male

Part #	Gas Service	Description
Z5004-4	Nitrous Oxide	QC/DISS x 1/4 NPT Male
Z5012-4	Nitrogen	QC/DISS x 1/4 NPT Male
Z5016-4	Air	QC/DISS x 1/4 NPT Male
Z5022-4	Vacuum	QC/DISS x 1/4 NPT Male
Not available	Oxygen	QC/DISS x 1/4 NPT Male
Z5032-4	Evacuation	QC/DISS x 1/4 NPT Male

DISS x 1/4 hose barb

Part #	Gas Service	Description
Z5104-4B	Nitrous Oxide	QC/DISS x 1/4 hose barb
Z5112-4B	Nitrogen	QC/DISS x 1/4 hose barb
Z5116-4B	Air	QC/DISS x 1/4 hose barb
Z5122-4B	Vacuum	QC/DISS x 1/4 hose barb
Z5124-4B	Oxygen	QC/DISS x 1/4 hose barb
Z5132-4B	Evacuation	QC/DISS x 1/4 hose barb
Z5122-5B	Vacuum	QC/DISS x 5/16 hose barb
Z5132-5B	Evacuation	QC/DISS x 5/16 hose barb

DISS Quick Connect Female x DISS Male

(used to adapt standard DISS nuts to Quick Connect)

Part #	Gas Service	Description
Z5204	Nitrous Oxide	QC/DISS x DISS male
Z5212	Nitrogen	QC/DISS x DISS male
Z5216	Air	QC/DISS x DISS male
Z5222	Vacuum	QC/DISS x DISS male
Z5224	Oxygen	QC/DISS x DISS male
Z5232	Evacuation	QC/DISS x DISS male